新冠抗原检测试剂产品手册
**Product Manual of
COVID-19 Antigen Detection Kit**

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COMPANY PROFILE

New Gene (Hangzhou) Bioengineering Co., Ltd. is located in Hangzhou, China. It is a high-tech company engaged in the research, development, manufacture and distribution of biological products. It is committed to creating biological materials such as antigens and antibodies, in vitro diagnostic reagents and related devices, and also the complete industrial chain of artificial intelligence assisted diagnosis system. The product line covers a full range of in vitro diagnostic products such as immune diagnosis, molecular diagnosis, and microbiological testing. NEWGENE has profound technical accumulation and unique technological advantages in the areas of early cancer screening, rapid detection of infectious diseases, and rapid screening of geriatric diseases.

NEWGENE's manufacturing system meets GMP standards for medical devices, and is certified with ISO13485 by British BSI. Relevant in vitro diagnostic reagent products have obtained the EU CE certification. NEWGENE is also a member on the "allow list" issued by Chinese Ministry of Commerce for anti-epidemic products exporting.

At present, NEWGENE COVID-19 Antigen Detection Kit has **registered in** many countries, including **Germany, France, Italy, Switzerland, Belgium, Portugal, Czech, Denmark, Hungary, Greece, Poland, Sweden, Moldova, Peru, Argentina, Ecuador, Kenya, Zimbabwe, Malaysia etc.**, and **Selftest of Germany, France, Czech, Denmark, Belgium, Sweden, Austria**, and passed the clinical **validation in** national lab in **Germany, Switzerland, Malaysia, Ecuador, Zimbabwe etc.** The products show good performance in sensitivity and specificity compared with international brand products and have exported to more than 50 countries and regions.

NEWGENE
Bioengineering

COVID-19 Antigen Detection Kit

Multiple Sampling Methods

Nasal Swab / Nasopharyngeal Swab / Oropharyngeal Swab / Sputum (Saliva)

Multiple Packaging Specifications

25 Tests/Box, 5 Tests/Box or 1 Test/Box

Global Recognition

Registered in 20+ Countries, Exported to 50+ Countries

Multiple Usage Scenarios

Professional Use & Home Use (Self-Testing)

Fast Detection

Results in 15 minutes

Superior Performance

High Sensitivity & Specificity

Contact Info

New Gene (Hangzhou) Bioengineering Co., Ltd.

Website: www.new-gene.com

www.new-gene.net

Email: marketing@new-gene.com

24-hour Hotline: (+86) 0571-5651 5020

COVID-19 Antigen Detection Kit (Nasal Swab Sample)

N0.	Components	25 Tests/Box	5 Tests/Box	1 Test/Box
1	Test Card	25	5	1
2	Sample Extraction Tube & Tube Cap	25	5	1
3	Sampling Swab: <i>for Nasal Swab</i>	25	5	1
4	Package Insert	1	1	1

25 Tests/Box

5 Tests/Box


1 Test/Box


Nasal Swab


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
97.1%	99.2%

COVID-19 Antigen Detection Kit (Nasopharyngeal Swab Sample)

NO.	Components	25 Tests/Box	1 Test/Box
1	Test Card	25	1
2	Sample Extraction Tube & Tube Cap	25	1
3	Sampling Swab: <i>for Nasopharyngeal Swab</i>	25	1
4	Package Insert	1	1

25 Tests/Box


1 Test/Box


Nasopharyngeal Swab


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
98.0%	99.1%

COVID-19 Antigen Detection Kit (Oropharyngeal Swab Sample)

N0.	Components	25 Tests/Box	1 Test/Box
1	Test Card	25	1
2	Sample Extraction Tube & Tube Cap	25	1
3	Sampling Swab: <i>for Oropharyngeal Swab</i>	25	1
4	Package Insert	1	1

25 Tests/Box

1 Test/Box

NEWGENE Bioengineering
COVID-19 Antigen Detection Kit
 Package Insert

CE IVD

Cat: COVID-19-NG28
 Version: E21-10-05
 Effective Date: 2021-05

System: Oropharyngeal Swab
 Effective Date: 2021-05

For professional use in vitro diagnostic use only.

PRODUCT NAME:
 COVID-19 Antigen Detection Kit

PACKING:
 1 box/each, 25 tests/box or 1 test/box.


INTENDED USE:
 This product is suitable for the qualitative detection of viral nucleic acid in oropharyngeal swab. It provides an aid in the diagnosis of infection with acute coronavirus.

SUMMARY:
 The novel coronavirus being the 21st genus COVID-19 is an acute respiratory infectious disease. It is generally asymptomatic. Currently, the patients infected by the novel coronavirus are the acute onset of infection, respiratory virus carriers can also be detected. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 5 days. The main transmission route are direct and indirect contact, respiratory, many cases are throat, sputum and diarrhea are also found in some cases.

PRINCIPLE:
 The COVID-19 Antigen Detection Kit is an immunochromatographic test kit that can rapidly detect nucleic acid in oropharyngeal swab. The test step is composed of the following steps: sample collection and storage, sample extraction, and detection. The reagent pad contains the colloidal gold conjugated with the monoclonal antibody against the nucleocapsid protein of SARS-CoV-2. The reaction membrane contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. The whole strip is fixed inside a plastic device. When the sample is added into the sample well, complete attachment or reaction will occur between the sample and the colloidal gold conjugate. If SARS-CoV-2 antigen is present in the sample, the complex of the anti-SARS-CoV-2 conjugate and the virus will be captured by the specific anti-SARS-CoV-2 monoclonal antibodies coated on the test line region. The absence of the test region's reaction result. To serve as a procedural control, a red line will always appear in the control line region (C) indicating that proper volume of sample has been added and membrane wetting effect has occurred.


COMPOSITION:
 1. Test Card
 2. Sample Extraction Tube
 3. Tube Cap
 4. Sampling Swab

STORAGE AND STABILITY:
 1. Store the product package at temperature 2-30°C or 38-89°F, and avoid exposure to sunlight. This kit is stable within the expiration date printed on the labeling.
 2. Once an aluminum foil pouch is opened, the test card inside should be used within one hour. Photocopy exposure to heat and humid environment may


Oropharyngeal Swab


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

Performance

Sensitivity	Specificity
95.7%	99.0%

COVID-19 Antigen Detection Kit (Sputum / Saliva Sample)

NO.	Components	25 Tests/Box	1 Test/Box
1	Test Card	25	1
2	Sample Extraction Tube & Tube Cap	25	1
3	Paper Cup	25	1
4	Sputum Dropper	25	1
5	Package Insert	1	1

25 Tests/Box


1 Test/Box


Sputum / Saliva


Test Procedures: Sampling

Test Procedures: Detection

Interpretation of Results

Positive(+): Red bands appear at both of T and C line in 15 to 30 minutes. A white band at the T line should be considered as a negative result.

Negative(-): A red band appears at C line while no red band appears at T line in 15 to 30 minutes after sample loading.

Invalid: As long as no red band appears at C line, it indicates that the test result is invalid, and should retest the sample with another test card.

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
A handwritten signature in blue ink, appearing to read 'M.J. van de Velde'.

Dr. M.J. van de Velde

DECLARATION OF CONFORMITY

Regarding In Vitro Diagnostic Directive (98/79/EC)

Manufacturer: New Gene (Hangzhou) Bioengineering Co., Ltd.

Address: Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou City, Zhejiang Province, P. R. China

EC Representative: SUNGO Europe B.V.

Address: Olympisch Stadion 24, 1076DE Amsterdam, Netherlands

Product Name: COVID-19 Antigen Detection Kit

Product Code: COVID-19-NG08

Specification: 25Tests/Box 1Test/Box

Classification: Others (IVDD)

Conformity Assessment

Procedure: Annex III of In Vitro Diagnostic Directive (98/79/EC)

We herewith declare that the above-mentioned products meet the requirements of In Vitro Diagnostic Directive (98/79/EC) and the following harmonized standards.

EN 23640:2015 EN 13640:2002

EN 980:2016 EN 13641:2002

EN ISO 14971:2019 EN ISO 18113-1:2011

EN 13612:2002 EN ISO 18113-4:2011

On behalf of SUNGO Europe office, I confirmed we are EU REP of the company who issue this document.

Signature:

Name/ Position: Mingfu Li / General Manager

Date: 29/09/2020

Place: Hangzhou, Zhejiang, China

Authorized Signature (S)

By Royal Charter

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2016

This is to certify that: New Gene (Hangzhou) Bioengineering Co., Ltd.
Room 1606, 16th Floor, No.5 Building
688 Bin'an Road
Binjiang District
Hangzhou
Zhejiang
310052
China

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Holds Certificate No: **MD 729179**

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Design and Development, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Drug Abuse, Manufacture and Distribution of In-vitro Diagnostic Rapid Test Kit of Infectious Diseases.

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For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

Original Registration Date: 2020-07-27

Latest Revision Date: 2020-07-27

Effective Date: 2020-07-27

Expiry Date: 2023-07-26

Page: 1 of 1

...making excellence a habit.™

Pass the evaluation by PEI of Germany

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Paul-Ehrlich-Institut 

26.03.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Aim

Comparison of different antigen rapid tests with using identical sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable CT values diluted 1:10 in negative samples in PBS. The CT values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standards. In the case of the PCRs used, a CT value of 25 corresponds to around 10^6 RNA copies/mL. 18 samples each were analysed with CT<25, 23 samples with CT between 25 and 30, and 9 samples with CT>30. The replication of the virus in cell culture was determined as a possible correlate for infectiousness as another characteristic of the samples.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 µL of the pool were analysed using the components of the test provided, e.g. swabs. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers with the same sample set allows an overview of the current state of art regarding sensitivity. The results do not allow any conclusions regarding specificity of the tests.

Those POCTs which have up to now been included in the evaluation and have been assessed as reflecting the current state of the art are listed in the table below. Other tests, which were assessed as not reflecting the state of the art were deleted from the list of the Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM). This comparative evaluation is constantly continued, and the table is amended accordingly.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for refunding, and that few other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

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Pass the evaluation by BAG of Switzerland

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Sars-CoV-2 Antigen Test

with

COVID-19 Antigen Detection Kit (Newgene(Hangzhou) Bioengineering)

Department of microbiology LMZ Dr Risch Pregassona

Validation interval:01/21 – 02/21

5 Conclusions

The comparable results with the R ***) Test makes the newgene test suitable to be introduced in Switzerland. It important to note the better sensitivity of that assay in comparison to R ***)

Pass the evaluation by MDA of Malaysia

INSTITUTE FOR MEDICAL RESEARCH KUALA LUMPUR

Performance of Newgene Bioengineering COVID-19 Antigen Detection Kit

Intended use

This product is suitable for the qualitative detection of novel coronavirus in nasopharyngeal swab or sputum samples. It provides an aid in the diagnosis of infection with novel coronavirus.

Manufacturer

New Gene (Hangzhou) Bioengineering Co., Ltd., China.

Test Principle

This kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect nucleocapsid from protein from SARS-CoV-2. When the sample is added into the sample well, colloidal-gold conjugated with the monoclonal antibody against the nucleocapsid protein of SARS-CoV-2 absorbed in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen is present in the sample, the complex of the anti-SARS-CoV-2 conjugate and the virus will be captured by the specific anti-SARS-CoV-2 monoclonal antibodies coated on the test line region (T). Absence of T line suggests a negative result. To serve as a procedural control, a red line will always appear in the control line region (C).

Test Kit

The evaluation was carried out using this kit with the lot number of 20210228-01, and the expiry date is Feb 27, 2022.

Instrument Used

NA

Reagent and Sample Preparation, Result Interpretation

Kindly refer to product package insert in the attachment

Sample Used

Sample type = Deep throat saliva (DTS)

Known Positive Sample= 50

Known Negative Sample= 50



Total samples used for analysis = 100

Performance Analysis

Test		Tested Kit Assay		Interpretation
		Positive	Negative	
IMR In-house Panel	SARS-CoV2 Positive (COVID-19)	48	2	Sensitivity= 96%
	SARS-CoV2 Negative (COVID-19)	0	50	Specificity = 100%

Comments

This kit was tested using DTS samples that had been tested positive using our RT-qPCR test system and the Ct values were ranged from 12 to 28. Saliva samples that had been tested negative using COVID-19 RT-qPCR test system and cultured samples of Influenza A H3, Influenza A pdm09, Adenovirus and Dengue were selected as negative panels.

Test performed by		Reviewed by	Approved by
			 3/5/22
Pn. E. Kavithambigai A/P Ellan Evaluator 1/Officer in charge	Pn. Hariyati Md. Ali Evaluator 2 /Officer in charge	Dr. Ravindran A/L Thayan Head of Virology Unit, IDRC, IMR	Dr. Hj. Tahir bin Aris Director of IMR

Disclaimer:

The content of this report does not imply that this product is endorsed or recommended by the Institute for Medical Research (IMR). This report may not be reviewed, abstracted, quoted, reproduced or translated, in part or in whole, without the prior written permission of IMR. IMR does not allow the use of this report or reference in any manner in the labeling or advertising of this product or kit package. However, the permission to import, export, or place in the market is subject to Medical Device Authority, Ministry of Health Malaysia's approval.

Pass the evaluation by Zimbabwe

NATIONAL MICROBIOLOGY REFERENCE LABORATORY

NEWGENECOVID 19 RAPID ANTIGEN TEST EVALUATION
AGAINST COVID 19 RT-PCR
COMPARISON REPORT

Conclusion

There was 100 % positive percent agreement between the NEWGENE SARS-CoV- 2 antigen test and SARS CoV-2 PCR test done using stored samples at NMRL. NEWGENE kit is recommended for use in Zimbabwe.

Evaluation team

Stanford Mupandasekwa	Lab Scientist	Signature <u>[Signature]</u>	Date <u>26/2/2021</u>
Boniface Muzividzi	Lab Scientist	Signature <u>[Signature]</u>	Date <u>26/2/2021</u>

Reviewed by:

Lucia Sisya	Quality officer	Signature <u>[Signature]</u>	Date <u>26/02/2021</u>
Agnes Juru	Chief lab Scientist	Signature <u>[Signature]</u>	Date <u>26/2/21</u>

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Registration or Allowed List (Partially)

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Self-Test of Denmark

New Gene (Hangzhou) Bioengineering Co., Ltd.
Building 5, 688 Bin'an Road, Changhe Subdistrict, Hangzhou, Zhejiang,
P.R. China

8 April 2021
Case no. 2021040180
T +45 44889775
E jast@dkma.dk

The Danish Medicines Agency hereby authorises New Gene (Hangzhou) Bioengineering to place on the market/put into service the COVID-19 Antigen Detection Kit in Denmark.

The authorisation is conditional with the following terms and conditions.

Terms and conditions:

- This authorisation is valid until 01.10.2021, or until the date when the device is **CE marked**, if this occurs before 01.10.2021. The manufacturer must inform the Danish Medicines Agency when the CE-mark has been acquired.
- The authorisation is limited to marketing/putting into service in Danish schools and Danish educational institutions.
- If the CE marking has not been achieved 01.07.2021 the manufacturer shall provide a status for the CE marking process.
- **New Gene Bioengineering** (manufacturer) shall establish a vigilance system, that ensures the collection and reporting to **The Danish Medicines Agency** of any incident, which occurs during the authorisation period.
- The manufacturer must ensure traceability of products through, for example, batch / LOT numbers.
- The information on the labelling and in the instruction for use (IFU) for the "COVID-19 Antigen Detection Kit" shall be in Danish, before distribution to the end-user.
- The labelling must clearly state that the products are placed on the market under this authorisation as a derogation from the conformity assessment procedures.
- The test must be accompanied by a detailed IFU, intended specifically for lay people and the intended user group (Danish educational institutions under supervision).
- Technical documentation for the product must be kept for a minimum of 5 years.
- A status shall be provided 01.07.2021 and by the end of this authorisation period of how many devices has been placed on the market under this authorisation.

The authorisation covers the following product:

Product: **COVID-19 Antigen Detection Kit**
Catalogue/ Model number(s): COVID-19-NG08
In Vitro Medical Device: General IVD, non-CE-marked for self-testing
Scope: immunoassay strip intended for the detection of SARS-CoV-2 nucleocapsid protein in human nasal swab samples.

Danish Medicines
Agency
Axel Heides Gade 1
2300 Copenhagen S
Denmark
T +45 44 88 95 95
E dkma@dkma.dk
LMST.DK

Get approved and listed by Belgium famhp as self-tests.

The screenshot shows the FAMHP website interface. At the top, there are language options (nl, fr, en) and a search bar. The main navigation menu includes 'Human use', 'Veterinary use', 'Information for the public', and 'Information for professionals'. The 'Human use' menu is currently selected. Below the navigation, there is a breadcrumb trail: Home > Human use > Health Products > Medical devices and their accessories > COVID-19 > Tests. The main content area is titled 'Tests' and contains several links: 'Antibody and antigen tests for professional use', 'New validation procedure for serological tests and antigenic tests', 'List of the recommended tests', 'Self-tests', 'List of SARS-CoV-2 antigen self-tests with a CE-certificate that can be sold in pharmacies to private persons', 'Procedure for manufacturers to market SARS-CoV-2 antigen rapid tests for professional use as self-tests in Belgium (see Dutch or French)', and 'List of professional rapid tests that may be sold as self-tests'. On the right side, there are two call-to-action boxes: 'Notification of adverse reactions or incidents' and 'PIL and SPC of a medicine'.

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Dnr: 5.2.3-2021-029480

New Gene Bioengineering.

Rm. 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street, Binjiang District, Hangzhou, Zhejiang, P. R. China

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GreMDIS - Greek Medical Devices Information System

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Product code, is the unique 13 digit code that a product acquires after its registration into the EOF's database

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MINISTRY OF HEALTH

PHARMACY AND POISONS BOARD

(Section 3B(2)(e) of the Pharmacy and Poisons Act, Cap 244 Laws of Kenya)

IN-VITRO DIAGNOSTIC EMERGENCY USE AUTHORIZATION

This Emergency Use Authorization is issued to **New Gene (Hangzhou) Bioengineering Co., Ltd.**, for distribution and sale of **Novel Coronavirus Antigen Detection Kit (Colloidal Gold)**

Emergency use Authorization (EUA) No.	MD/2021/7674
EUA valid until	End of COVID -19 Pandemic or EUA revocation
Device category	Medical Device class C/D
GMDN	N/A
GMDN Term	N/A
Intended purpose	For epidemiological COVID-19 Screening
Conditional Approval	N/A

EUA No.: MD/2021/7674

Date of Authorization: March 4, 2021

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Razón Social: **ALCAT S.A.**